

For Research Use Only

HSD17B4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51363-1

Capture Antibody Information

Catalog Number:
60951-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1A3H11

Reactivity:
human

GenBank:
BC003098

Immunogen Catalog Number:
Ag7221

Conjugate:
Unconjugated

Full name:
hydroxysteroid (17-beta) dehydrogenase 4

Gene ID:
3295

Detection Antibody Information

Catalog Number:
60951-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B11F1

Reactivity:
human

GenBank:
BC003098

Immunogen Catalog Number:
Ag7221

Conjugate:
Unconjugated

Full name:
hydroxysteroid (17-beta) dehydrogenase 4

Gene ID:
3295

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51363-1 targets HSD17B4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSD17B4 Monoclonal antibody, PBS Only (Capture) 60951-1-PBS (1A3H11). 100 µg. Concentration 1 mg/mL.

Detection antibody: HSD17B4 Monoclonal antibody, PBS Only (Detector) 60951-2-PBS (1B11F1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

