

For Research Use Only

GFM2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51366-1

Capture Antibody Information

Catalog Number:
60952-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G1F8

Reactivity:
human

GenBank:
BC015712

Immunogen Catalog Number:
Ag36976

Conjugate:
Unconjugated

Full name:
G elongation factor, mitochondrial 2

Gene ID:
84340

Detection Antibody Information

Catalog Number:
60952-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H1G1

Reactivity:
human

GenBank:
BC015712

Immunogen Catalog Number:
Ag36976

Conjugate:
Unconjugated

Full name:
G elongation factor, mitochondrial 2

Gene ID:
84340

Applications

Tested Applications:
Cytometric bead array

Range:
12.5-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51366-1 targets GFM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GFM2 Monoclonal antibody, PBS Only (Capture) 60952-1-PBS (1G1F8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: GFM2 Monoclonal antibody, PBS Only (Detector) 60952-2-PBS (1H1G1). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51366-1, GFM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60952-1-PBS. Detection antibody: 60952-2-PBS.



Cytometric bead array standard curve of MP51366-1, GFM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60952-1-PBS. Detection antibody: 60952-2-PBS. Standard: Ag36976. Range: 12.5-400 ng/mL.