

For Research Use Only

TAF11 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51378-1

Capture Antibody Information

Catalog Number:
60962-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3B1G7
Reactivity:
human
GenBank:
BC021972
Immunogen Catalog Number:
Ag9629

Conjugate:
Unconjugated
Full name:
TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa
Gene ID:
6882

Detection Antibody Information

Catalog Number:
60962-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
3B10G9
Reactivity:
human
GenBank:
BC021972
Immunogen Catalog Number:
Ag9629

Conjugate:
Unconjugated
Full name:
TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa
Gene ID:
6882

Applications

Tested Applications:
Cytometric bead array

Range:
12.5-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51378-1 targets TAF11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAF11 Monoclonal antibody, PBS Only (Capture) 60962-1-PBS (3B1G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: TAF11 Monoclonal antibody, PBS Only (Detector) 60962-2-PBS (3B10G9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS.

Cytometric bead array standard curve of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS. Standard: Ag9629. Range: 12.5-200 ng/mL.