

DTX2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51381-1**

Capture Antibody Information

Catalog Number:
67209-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
2D4F8
Reactivity:
human
GenBank:
BC064535
Immunogen Catalog Number:
Ag12022

Conjugate:
Unconjugated
Full name:
deltex homolog 2 (Drosophila)
Gene ID:
113878

Detection Antibody Information

Catalog Number:
67209-2-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
1B8E2
Reactivity:
human
GenBank:
BC064535
Immunogen Catalog Number:
Ag12022

Conjugate:
Unconjugated
Full name:
deltex homolog 2 (Drosophila)
Gene ID:
113878

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51381-1 targets DTX2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DTX2 Monoclonal antibody, PBS Only (Capture) 67209-1-PBS (2D4F8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: DTX2 Monoclonal antibody, PBS Only (Detector) 67209-2-PBS (1B8E2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51381-1, DTX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67209-1-PBS. Detection antibody: 67209-2-PBS.



Cytometric bead array standard curve of MP51381-1, DTX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67209-1-PBS. Detection antibody: 67209-2-PBS. Standard:Ag12022. Range: 3.125-400 ng/mL