

FHL3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP51383-1

Capture Antibody Information

Catalog Number: 60966-1-PBS	Clone ID: 1G3A10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: four and a half LIM domains 3
Isotype: IgG2b	GenBank: BC011697	Gene ID: 2275
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag1495	

Detection Antibody Information

Catalog Number: 60966-2-PBS	Clone ID: 3D7A12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: four and a half LIM domains 3
Isotype: IgG2a	GenBank: BC011697	Gene ID: 2275
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag1495	

Applications

Tested Applications: Cytometric bead array	Range: 6.25-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP51383-1 targets FHL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FHL3 Monoclonal antibody, PBS Only (Capture) 60966-1-PBS (1G3A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: FHL3 Monoclonal antibody, PBS Only (Detector) 60966-2-PBS (3D7A12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS.



Cytometric bead array standard curve of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS. Standard: Ag1495. Range: 6.25-400 ng/mL.