

For Research Use Only

STK4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51388-1

Capture Antibody Information

Catalog Number:
66663-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C4B11
Reactivity:
human
GenBank:
BC093768
Immunogen Catalog Number:
Ag17738

Conjugate:
Unconjugated
Full name:
serine/threonine kinase 4
Gene ID:
6789

Detection Antibody Information

Catalog Number:
66663-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3D1H11
Reactivity:
human
GenBank:
BC093768
Immunogen Catalog Number:
Ag17738

Conjugate:
Unconjugated
Full name:
serine/threonine kinase 4
Gene ID:
6789

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51388-1 targets STK4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STK4 Monoclonal antibody, PBS Only (Capture) 66663-2-PBS (2C4B11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: STK4 Monoclonal antibody, PBS Only (Detector) 66663-3-PBS (3D1H11). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

