

For Research Use Only

YPEL3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51392-1

Capture Antibody Information

Catalog Number: 60970-1-PBS	Clone ID: 1B7B11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: yippee-like 3 (Drosophila)
Isotype: IgG1	GenBank: BC005009	Gene ID: 83719
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7626	

Detection Antibody Information

Catalog Number: 60970-2-PBS	Clone ID: 1B5F4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: yippee-like 3 (Drosophila)
Isotype: IgG1	GenBank: BC005009	Gene ID: 83719
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag7626	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51392-1 targets YPEL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: YPEL3 Monoclonal antibody, PBS Only (Capture) 60970-1-PBS (1B7B11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: YPEL3 Monoclonal antibody, PBS Only (Detector) 60970-2-PBS (1B5F4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51392-1, YPEL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60970-1-PBS. Detection antibody: 60970-2-PBS.



Cytometric bead array standard curve of MP51392-1, YPEL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60970-1-PBS. Detection antibody: 60970-2-PBS. Standard: Ag7626. Range: 0.781-100 ng/mL