

For Research Use Only

LMNB2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51397-1

Capture Antibody Information

Catalog Number:
68475-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2B6H1

Reactivity:
human

GenBank:
BC006551

Immunogen Catalog Number:
Ag33296

Conjugate:
Unconjugated

Full name:
lamin B2

Gene ID:
84823

Detection Antibody Information

Catalog Number:
68475-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E9A5

Reactivity:
human

GenBank:
BC006551

Immunogen Catalog Number:
Ag33296

Conjugate:
Unconjugated

Full name:
lamin B2

Gene ID:
84823

Applications

Tested Applications:
Cytometric bead array

Range:
0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51397-1 targets LMNB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LMNB2 Monoclonal antibody, PBS Only (Capture) 68475-2-PBS (2B6H1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LMNB2 Monoclonal antibody, PBS Only (Detector) 68475-3-PBS (1E9A5). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51397-1, LMNB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68475-2-PBS. Detection antibody: 68475-3-PBS. Standard: Ag33296. Range: 0.195-100 ng/mL.

Sample test of MP51397-1, LMNB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68475-2-PBS. Detection antibody: 68475-3-PBS.