

For Research Use Only

# SYNJ1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51398-1

## Capture Antibody Information

**Catalog Number:**  
60973-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1B5C5  
**Reactivity:**  
human  
**GenBank:**  
BC136603  
**Immunogen Catalog Number:**  
Ag17719

**Conjugate:**  
Unconjugated  
**Full name:**  
synaptojanin 1  
**Gene ID:**  
8867

## Detection Antibody Information

**Catalog Number:**  
60973-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1C8B1  
**Reactivity:**  
human  
**GenBank:**  
BC136603  
**Immunogen Catalog Number:**  
Ag17719

**Conjugate:**  
Unconjugated  
**Full name:**  
synaptojanin 1  
**Gene ID:**  
8867

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51398-1 targets SYNJ1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: SYNJ1 Monoclonal antibody, PBS Only (Capture) 60973-1-PBS (1B5C5). 100 µg. Concentration 1 mg/mL.  
Detection antibody: SYNJ1 Monoclonal antibody, PBS Only (Detector) 60973-2-PBS (1C8B1). 100 µg. Concentration 1 mg/mL.  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.



## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sample test of MP51398-1, SYNJ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60973-1-PBS. Detection antibody: 60973-2-PBS.

Cytometric bead array standard curve of MP51398-1, SYNJ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60973-1-PBS. Detection antibody: 60973-2-PBS. Standard: Ag17719. Range: 0.195-50 ng/mL.