

TCEB3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51401-1

Capture Antibody Information

Catalog Number:
60974-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B1A9
Reactivity:
human
GenBank:
BC002883
Immunogen Catalog Number:
Ag33045

Conjugate:
Unconjugated
Full name:
transcription elongation factor B (SIII), polypeptide 3 (110kDa, elongin A)
Gene ID:
6924

Detection Antibody Information

Catalog Number:
60974-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2B1A3
Reactivity:
human
GenBank:
BC002883
Immunogen Catalog Number:
Ag33045

Conjugate:
Unconjugated
Full name:
transcription elongation factor B (SIII), polypeptide 3 (110kDa, elongin A)
Gene ID:
6924

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51401-1 targets TCEB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TCEB3 Monoclonal antibody, PBS Only (Capture) 60974-1-PBS (1B1A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: TCEB3 Monoclonal antibody, PBS Only (Detector) 60974-2-PBS (2B1A3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51401-1, TCEB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60974-1-PBS. Detection antibody: 60974-2-PBS.

Cytometric bead array standard curve of MP51401-1, TCEB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60974-1-PBS. Detection antibody: 60974-2-PBS. Standard: Ag33045. Range: 0.098-100 ng/mL