

CRX Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51402-2**

Capture Antibody Information

Catalog Number: 67555-1-PBS	Clone ID: 1C3C11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat, pig	Full name: cone-rod homeobox
Isotype: IgG1	GenBank: BC016664	Gene ID: 1406
Purification Method: Protein G purification	Immunogen Catalog Number: Ag30086	

Detection Antibody Information

Catalog Number: 67555-4-PBS	Clone ID: 2H3C9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cone-rod homeobox
Isotype: IgG2b	GenBank: BC016664	Gene ID: 1406
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag30086	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51402-2 targets CRX in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CRX Monoclonal antibody, PBS Only (Capture) 67555-1-PBS (1C3C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CRX Monoclonal antibody, PBS Only (Detector) 67555-4-PBS (2H3C9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51402-2, CRX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67555-1-PBS. Detection antibody: 67555-4-PBS.



Cytometric bead array standard curve of MP51402-2, CRX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67555-1-PBS. Detection antibody: 67555-4-PBS. Standard: Ag30086. Range: 0.098-100 ng/mL