

# VSTM2L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51415-1

## Capture Antibody Information

**Catalog Number:**  
60981-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1D11H9  
**Reactivity:**  
human  
**GenBank:**  
BC033818  
**Immunogen Catalog Number:**  
Ag22362

**Conjugate:**  
Unconjugated  
**Full name:**  
V-set and transmembrane domain containing 2 like  
**Gene ID:**  
128434

## Detection Antibody Information

**Catalog Number:**  
60981-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1B6H3  
**Reactivity:**  
human  
**GenBank:**  
BC033818  
**Immunogen Catalog Number:**  
Ag22362

**Conjugate:**  
Unconjugated  
**Full name:**  
V-set and transmembrane domain containing 2 like  
**Gene ID:**  
128434

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
6.25-400 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51415-1 targets VSTM2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VSTM2L Monoclonal antibody, PBS Only (Capture) 60981-1-PBS (1D11H9). 100 µg. Concentration 1 mg/mL.

Detection antibody: VSTM2L Monoclonal antibody, PBS Only (Detector) 60981-2-PBS (1B6H3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

