

For Research Use Only

TAF13 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51418-1

Capture Antibody Information

Catalog Number:
60984-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2G4C12

Reactivity:
human

GenBank:
BC121180

Immunogen Catalog Number:
Ag18502

Conjugate:
Unconjugated

Full name:
TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa

Gene ID:
6884

Detection Antibody Information

Catalog Number:
60984-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2A1D8

Reactivity:
human

GenBank:
BC121180

Immunogen Catalog Number:
Ag18502

Conjugate:
Unconjugated

Full name:
TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa

Gene ID:
6884

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51418-1 targets TAF13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAF13 Monoclonal antibody, PBS Only (Capture) 60984-1-PBS (2G4C12). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TAF13 Monoclonal antibody, PBS Only (Detector) 60984-2-PBS (2A1D8). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51418-1, TAF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60984-1-PBS. Detection antibody: 60984-2-PBS. Standard: Ag18502. Range: 0.098-100 ng/mL.

Cytometric bead array sample test of MP51418-1, TAF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60984-1-PBS. Detection antibody: 60984-2-PBS.