

SERPINA12 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51421-1**

Capture Antibody Information

Catalog Number: 67573-1-PBS	Clone ID: 2A1B6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 12
Isotype: IgG1	GenBank: BC040857	Gene ID: 145264
Purification Method: Protein A purification	Immunogen Catalog Number: Ag30136	

Detection Antibody Information

Catalog Number: 67573-2-PBS	Clone ID: 1F9C5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 12
Isotype: IgG1	GenBank: BC040857	Gene ID: 145264
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag30136	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51421-1 targets SERPINA12 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Serpin A12/Vaspin Monoclonal antibody, PBS Only (Capture) 67573-1-PBS (2A1B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: SERPINA12 Monoclonal antibody, PBS Only (Detector) 67573-2-PBS (1F9C5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

