

CAPRIN1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51429-1**

Capture Antibody Information

Catalog Number: 66352-1-PBS	Clone ID: 6B10F11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat	Full name: cell cycle associated protein 1
Isotype: IgG1	GenBank: BC001731	Gene ID: 4076
Purification Method: Protein G purification	Immunogen Catalog Number: Ag7413	

Detection Antibody Information

Catalog Number: 66352-2-PBS	Clone ID: 2C7B2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cell cycle associated protein 1
Isotype: IgG2a	GenBank: BC001731	Gene ID: 4076
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag7413	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51429-1 targets CAPRIN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CAPRIN1 Monoclonal antibody, PBS Only (Capture) 66352-1-PBS (6B10F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CAPRIN1 Monoclonal antibody, PBS Only (Detector) 66352-2-PBS (2C7B2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51429-1, CAPRN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66352-1-PBS. Detection antibody: 66352-2-PBS. Standard:Ag7413. Range: 1.563-200 ng/mL.



Cytometric bead array sample test of MP51429-1, CAPRN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66352-1-PBS. Detection antibody: 66352-2-PBS.