

TAU Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51432-1**

Capture Antibody Information

Catalog Number:
60995-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
4D7F3
Reactivity:
human
GenBank:
BC000558

Conjugate:
Unconjugated
Full name:
microtubule-associated protein tau
Gene ID:
4137

Detection Antibody Information

Catalog Number:
66499-2-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A Magarose purification

Clone ID:
1D7C9
Reactivity:
human
GenBank:
BC000558
Immunogen Catalog Number:
Ag21926

Conjugate:
Unconjugated
Full name:
microtubule-associated protein tau
Gene ID:
4137

Applications

Tested Applications:
Cytometric bead array

Range:
Physical level (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51432-1 targets TAU in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAU Monoclonal antibody, PBS Only (Capture) 60995-1-PBS (4D7F3). 100 µg. Concentration 1 mg/mL.

Detection antibody: TAU Monoclonal antibody, PBS Only (Detector) 66499-2-PBS (1D7C9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51432-1, TAU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60995-1-PBS. Detection antibody: 66499-2-PBS.