

For Research Use Only

# DHX15 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51437-1

## Capture Antibody Information

Catalog Number:  
60999-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1F7C5

Reactivity:  
human

GenBank:  
BC035974

Immunogen Catalog Number:  
Ag23433

Conjugate:  
Unconjugated

Full name:  
DEAH (Asp-Glu-Ala-His) box polypeptide 15

Gene ID:  
1665

## Detection Antibody Information

Catalog Number:  
60999-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2A2F8

Reactivity:  
human

GenBank:  
BC035974

Immunogen Catalog Number:  
Ag23433

Conjugate:  
Unconjugated

Full name:  
DEAH (Asp-Glu-Ala-His) box polypeptide 15

Gene ID:  
1665

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51437-1 targets DHX15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DHX15 Monoclonal antibody, PBS Only (Capture) 60999-1-PBS (1F7C5). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: DHX15 Monoclonal antibody, PBS Only (Detector) 60999-2-PBS (2A2F8). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51437-1, DHX15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60999-1-PBS. Detection antibody: 60999-2-PBS. Standard: Ag23433. Range: 0.391-200 ng/mL.

Cytometric bead array sample test of MP51437-1, DHX15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60999-1-PBS. Detection antibody: 60999-2-PBS.