

LIMA1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51439-1

Capture Antibody Information

Catalog Number:
61000-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3F12C6

Reactivity:
human

GenBank:
BC001247

Immunogen Catalog Number:
Ag13288

Conjugate:
Unconjugated

Full name:
LIM domain and actin binding 1

Gene ID:
51474

Detection Antibody Information

Catalog Number:
61000-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
4C8C7

Reactivity:
human

GenBank:
BC001247

Immunogen Catalog Number:
Ag13288

Conjugate:
Unconjugated

Full name:
LIM domain and actin binding 1

Gene ID:
51474

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51439-1 targets LIMA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LIMA1 Monoclonal antibody, PBS Only (Capture) 61000-1-PBS (3F12C6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LIMA1 Monoclonal antibody, PBS Only (Detector) 61000-2-PBS (4C8C7). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C


Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51439-1, LIMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61000-1-PBS. Detection antibody: 61000-2-PBS. Standard: Ag13288. Range: 1.563-200 ng/mL.



Cytometric bead array sample test of MP51439-1, LIMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61000-1-PBS. Detection antibody: 61000-2-PBS.