

RBM15 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51442-1**

Capture Antibody Information

Catalog Number:
66059-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A purification

Clone ID:
4A1A4
Reactivity:
human
GenBank:
BC006397
Immunogen Catalog Number:
Ag17970

Conjugate:
Unconjugated
Full name:
RNA binding motif protein 15
Gene ID:
64783

Detection Antibody Information

Catalog Number:
66059-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3H2B6
Reactivity:
human
GenBank:
BC006397
Immunogen Catalog Number:
Ag17970

Conjugate:
Unconjugated
Full name:
RNA binding motif protein 15
Gene ID:
64783

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51442-1 targets RBM15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBM15 Monoclonal antibody, PBS Only (Capture) 66059-1-PBS (4A1A4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RBM15 Monoclonal antibody, PBS Only (Detector) 66059-2-PBS (3H2B6). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51442-1, RBM15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66059-1-PBS. Detection antibody: 66059-2-PBS. Standard: Ag17970. Range: 0.391-100 ng/mL



Cytometric bead array sample test of MP51442-1, RBM15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66059-1-PBS. Detection antibody: 66059-2-PBS.