

BMI1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51443-1**

Capture Antibody Information

Catalog Number: 66161-2-PBS	Clone ID: 2E3E9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: BMI1 polycomb ring finger oncogene
Isotype: IgG2b	GenBank: BC011652	Gene ID: 648
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag21284	

Detection Antibody Information

Catalog Number: 66161-1-PBS	Clone ID: 2E10A11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, monkey	Full name: BMI1 polycomb ring finger oncogene
Isotype: IgG2b	GenBank: BC011652	Gene ID: 648
Purification Method: Protein A purification	Immunogen Catalog Number: Ag21284	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51443-1 targets BMI1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BMI1 Monoclonal antibody, PBS Only (Capture) 66161-2-PBS (2E3E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: BMI1 Monoclonal antibody, PBS Only (Detector) 66161-1-PBS (2E10A11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51443-1, BMI1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66161-2-PBS. Detection antibody: 66161-1-PBS. Standard:Ag21284. Range: 0.391-100 ng/mL.



Cytometric bead array sample test of MP51443-1, BMI1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66161-2-PBS. Detection antibody: 66161-1-PBS.