

For Research Use Only

CMIP Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51453-1

Capture Antibody Information

Catalog Number: 61004-1-PBS	Clone ID: 2A7A11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: c-Maf-inducing protein
Isotype: IgG1	GenBank: BC038113	Gene ID: 80790
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17875	

Detection Antibody Information

Catalog Number: 61004-2-PBS	Clone ID: 2B12A9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: c-Maf-inducing protein
Isotype: IgG1	GenBank: BC038113	Gene ID: 80790
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17875	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51453-1 targets CMIP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CMIP Monoclonal antibody, PBS Only (Capture) 61004-1-PBS (2A7A11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CMIP Monoclonal antibody, PBS Only (Detector) 61004-2-PBS (2B12A9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage



Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51453-1, CMIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61004-1-PBS. Detection antibody: 61004-2-PBS. Standard: Ag17875. Range: 1.563-400 ng/mL.

Cytometric bead array sample test of MP51453-1, CMIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61004-1-PBS. Detection antibody: 61004-2-PBS.