

CBLL1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51461-1**

Capture Antibody Information

Catalog Number: 61007-1-PBS	Clone ID: 1A10A2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1
Isotype: IgG1	GenBank: BC027460	Gene ID: 79872
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag30913	

Detection Antibody Information

Catalog Number: 61007-2-PBS	Clone ID: 2A1F7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1
Isotype: IgG1	GenBank: BC027460	Gene ID: 79872
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag30913	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51461-1 targets CBLL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CBLL1 Monoclonal antibody, PBS Only (Capture) 61007-1-PBS (1A10A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: CBLL1 Monoclonal antibody, PBS Only (Detector) 61007-2-PBS (2A1F7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51461-1, CBLL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61007-1-PBS. Detection antibody: 61007-2-PBS. Standard:Ag30913. Range: 3.125-400 ng/mL.



Cytometric bead array sample test of MP51461-1, CBLL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61007-1-PBS. Detection antibody: 61007-2-PBS.