

For Research Use Only

BTN3A1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51463-2

Capture Antibody Information

Catalog Number:
61008-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E6H10

Reactivity:
human

GenBank:
BC121800

Immunogen Catalog Number:
Eg1915

Conjugate:
Unconjugated

Full name:
butyrophilin, subfamily 3, member A1

Gene ID:
11119

Detection Antibody Information

Catalog Number:
61008-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D10E5

Reactivity:
human

GenBank:
BC121800

Immunogen Catalog Number:
Eg1915

Conjugate:
Unconjugated

Full name:
butyrophilin, subfamily 3, member A1

Gene ID:
11119

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51463-2 targets BTN3A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BTN3A1 Monoclonal antibody, PBS Only (Capture) 61008-1-PBS (1E6H10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: BTN3A1 Monoclonal antibody, PBS Only (Detector) 61008-3-PBS (2D10E5). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51463-2, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-3-PBS. Standard: Eg1915. Range: 0.391-50 ng/mL



Cytometric bead array sample test of MP51463-2, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-3-PBS. Standard: Eg1915. Range: 0.391-50 ng/mL