

CCDC6 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP51465-1

Capture Antibody Information

Catalog Number:
67637-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D6A8
Reactivity:
human, mouse, rat
GenBank:
BC036757
Immunogen Catalog Number:
Ag6952

Conjugate:
Unconjugated
Full name:
coiled-coil domain containing 6
Gene ID:
8030

Detection Antibody Information

Catalog Number:
67637-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F10F2
Reactivity:
human
GenBank:
BC036757
Immunogen Catalog Number:
Ag6952

Conjugate:
Unconjugated
Full name:
coiled-coil domain containing 6
Gene ID:
8030

Applications

Tested Applications:
Cytometric bead array

Range:
0.078-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51465-1 targets CCDC6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC6 Monoclonal antibody, PBS Only (Capture) 67637-1-PBS (1D6A8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CCDC6 Monoclonal antibody, PBS Only (Detector) 67637-2-PBS (2F10F2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51465-1, CCDC6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67637-1-PBS. Detection antibody: 67637-2-PBS.



Cytometric bead array standard curve of MP51465-1, CCDC6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67637-1-PBS. Detection antibody: 67637-2-PBS. Standard:Ag6952. Range: 0.078-10 ng/mL