

For Research Use Only

HSPG2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51466-1

Capture Antibody Information

Catalog Number:
61009-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F11B12
Reactivity:
human
GenBank:
NM_001291860
Immunogen Catalog Number:
Ag30846

Conjugate:
Unconjugated
Full name:
heparan sulfate proteoglycan 2
Gene ID:
3339

Detection Antibody Information

Catalog Number:
61009-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F6D5
Reactivity:
human
GenBank:
NM_001291860
Immunogen Catalog Number:
Ag30846

Conjugate:
Unconjugated
Full name:
heparan sulfate proteoglycan 2
Gene ID:
3339

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51466-1 targets HSPG2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HSPG2 Monoclonal antibody, PBS Only (Capture) 61009-1-PBS (1F11B12). 100 µg. Concentration 1 mg/mL.

Detection antibody: HSPG2 Monoclonal antibody, PBS Only (Detector) 61009-2-PBS (1F6D5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51466-1, HSPG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61009-1-PBS. Detection antibody: 61009-2-PBS.



Cytometric bead array standard curve of MP51466-1, HSPG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61009-1-PBS. Detection antibody: 61009-2-PBS. Standard: Ag30846. Range: 1.563-200 ng/mL