

OTUD7B Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51475-1**

Capture Antibody Information

Catalog Number: 66276-1-PBS	Clone ID: 2C12D9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: OTU domain containing 7B
Isotype: IgG2b	GenBank: BC020622	Gene ID: 56957
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10513	

Detection Antibody Information

Catalog Number: 66276-2-PBS	Clone ID: 1H6C5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: OTU domain containing 7B
Isotype: IgG1	GenBank: BC020622	Gene ID: 56957
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag10513	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51475-1 targets OTUD7B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OTUD7B Monoclonal antibody, PBS Only (Capture) 66276-1-PBS (2C12D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: OTUD7B Monoclonal antibody, PBS Only (Detector) 66276-2-PBS (1H6C5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51475-1, OTUD7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66276-1-PBS. Detection antibody: 66276-2-PBS.



Cytometric bead array standard curve of MP51475-1, OTUD7B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66276-1-PBS. Detection antibody: 66276-2-PBS. Standard:Ag10513. Range: 0.781-100 ng/mL