

EPHX2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51477-1**

Capture Antibody Information

Catalog Number: 60083-2-PBS	Clone ID: 5B8C12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: epoxide hydrolase 2, cytoplasmic
Isotype: IgG1	GenBank: BC013874	Gene ID: 2053
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag4870	

Detection Antibody Information

Catalog Number: 60083-1-PBS	Clone ID: 3A7C4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: epoxide hydrolase 2, cytoplasmic
Isotype: IgG1	GenBank: BC013874	Gene ID: 2053
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4870	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51477-1 targets EPX2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPX2 Monoclonal antibody, PBS Only (Capture) 60083-2-PBS (5B8C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: EPX2 Monoclonal antibody, PBS Only (Detector) 60083-1-PBS (3A7C4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51477-1, EPHX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60083-2-PBS. Detection antibody: 60083-1-PBS.



Cytometric bead array standard curve of MP51477-1, EPHX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60083-2-PBS. Detection antibody: 60083-1-PBS. Standard: Ag1283. Range: 1.563-400 ng/mL.