

ALPL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51482-1**

Capture Antibody Information

Catalog Number:
68719-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F5G3
Reactivity:
human
GenBank:
BC021289
Immunogen Catalog Number:
Ag33678

Conjugate:
Unconjugated
Full name:
alkaline phosphatase,
liver/bone/kidney
Gene ID:
249

Detection Antibody Information

Catalog Number:
68719-3-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
2C10H7
Reactivity:
human
GenBank:
BC021289
Immunogen Catalog Number:
Ag33678

Conjugate:
Unconjugated
Full name:
alkaline phosphatase,
liver/bone/kidney
Gene ID:
249

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51482-1 targets ALPL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ALPL Monoclonal antibody, PBS Only (Capture) 68719-2-PBS (1F5G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: ALPL Monoclonal antibody, PBS Only (Detector) 68719-3-PBS (2C10H7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51482-1, ALPL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68719-2-PBS. Detection antibody: 68719-3-PBS.



Cytometric bead array standard curve of MP51482-1, ALPL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68719-2-PBS. Detection antibody: 68719-3-PBS. Standard: Ag33678. Range: 0.781-100 ng/mL.