

CCAR2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51483-1**

Capture Antibody Information

Catalog Number: 66497-2-PBS	Clone ID: 2A10A7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: KIAA1967
Isotype: IgG1	GenBank: BC018269	Gene ID: 57805
Purification Method: Protein G purification	Immunogen Catalog Number: Ag18374	

Detection Antibody Information

Catalog Number: 66497-1-PBS	Clone ID: 2D10F8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, pig	Full name: KIAA1967
Isotype: IgG2a	GenBank: BC018269	Gene ID: 57805
Purification Method: Protein A purification	Immunogen Catalog Number: Ag18374	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51483-1 targets CCAR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DBC1 Monoclonal antibody, PBS Only (Capture) 66497-2-PBS (2A10A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCAR2 Monoclonal antibody, PBS Only (Detector) 66497-1-PBS (2D10F8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51483-1, CCAR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66497-2-PBS. Detection antibody: 66497-1-PBS.



Cytometric bead array standard curve of MP51483-1, CCAR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66497-2-PBS. Detection antibody: 66497-1-PBS. Standard: Ag18374. Range: 1.563-100 ng/mL.