

RABEPK Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51484-1**

Capture Antibody Information

Catalog Number:
66622-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1E11A3
Reactivity:
human, mouse, rat
GenBank:
BC065725
Immunogen Catalog Number:
Ag7796

Conjugate:
Unconjugated
Full name:
Rab9 effector protein with kelch motifs
Gene ID:
10244

Detection Antibody Information

Catalog Number:
66622-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F11A3
Reactivity:
human
GenBank:
BC065725
Immunogen Catalog Number:
Ag7796

Conjugate:
Unconjugated
Full name:
Rab9 effector protein with kelch motifs
Gene ID:
10244

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51484-1 targets RABEPK in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RABEPK/p40 Monoclonal antibody, PBS Only (Capture) 66622-1-PBS (1E11A3). 100 µg. Concentration 1 mg/mL.

Detection antibody: RABEPK Monoclonal antibody, PBS Only (Detector) 66622-2-PBS (1F11A3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51484-1, RABEPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66622-1-PBS. Detection antibody: 66622-2-PBS.



Cytometric bead array standard curve of MP51484-1, RABEPK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66622-1-PBS. Detection antibody: 66622-2-PBS. Standard:Ag2593. Range: 1.563-400 ng/mL.