

IKBKG Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51489-1**

Capture Antibody Information

Catalog Number:
66460-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F2E2
Reactivity:
human
GenBank:
BC012114
Immunogen Catalog Number:
Ag13358

Conjugate:
Unconjugated
Full name:
inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma
Gene ID:
8517

Detection Antibody Information

Catalog Number:
66460-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1A1A2
Reactivity:
human
GenBank:
BC012114
Immunogen Catalog Number:
Ag13358

Conjugate:
Unconjugated
Full name:
inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma
Gene ID:
8517

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51489-1 targets IKBKG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IKBKG Monoclonal antibody, PBS Only (Capture) 66460-1-PBS (1F2E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: IKBKG Monoclonal antibody, PBS Only (Detector) 66460-2-PBS (1A1A2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Selected Validation Data

