

RHOC Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51492-1**

Capture Antibody Information

Catalog Number:
67542-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2G4C8
Reactivity:
human
GenBank:
BC007245
Immunogen Catalog Number:
Ag14633

Conjugate:
Unconjugated
Full name:
ras homolog gene family, member C
Gene ID:
389

Detection Antibody Information

Catalog Number:
67542-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3E11B4
Reactivity:
human
GenBank:
BC007245
Immunogen Catalog Number:
Ag14633

Conjugate:
Unconjugated
Full name:
ras homolog gene family, member C
Gene ID:
389

Applications

Tested Applications:
Cytometric bead array

Range:
6.250-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51492-1 targets RHOC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RHOC Monoclonal antibody, PBS Only (Capture) 67542-2-PBS (2G4C8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RHOC Monoclonal antibody, PBS Only (Detector) 67542-3-PBS (3E11B4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51492-1, RHOC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67542-2-PBS. Detection antibody: 67542-3-PBS.



Cytometric bead array standard curve of MP51492-1, RHOC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67542-2-PBS. Detection antibody: 67542-3-PBS. Standard:Ag14633. Range: 6.250-400 ng/mL.