

CCNE2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51494-1**

Capture Antibody Information

Catalog Number: 61015-1-PBS	Clone ID: 1A7H7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cyclin E2
Isotype: IgG1	GenBank: BC020729	Gene ID: 9134
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag23601	

Detection Antibody Information

Catalog Number: 61016-1-PBS	Clone ID: 1E7B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: cyclin E2
Isotype: IgG3	GenBank: BC020729	Gene ID: 9134
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29625	

Applications

Tested Applications: Cytometric bead array	Range: 6.250-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51494-1 targets CCNE2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Cyclin E2 Monoclonal antibody, PBS Only (Capture) 61015-1-PBS (1A7H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Cyclin E2 Monoclonal antibody, PBS Only (Detector) 61016-1-PBS (1E7B5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51494-1, CCNE2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61015-1-PBS. Detection antibody: 61016-1-PBS.



Cytometric bead array standard curve of MP51494-1, CCNE2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61015-1-PBS. Detection antibody: 61016-1-PBS. Standard:Ag23601. Range: 6.250-200 ng/mL.