

U2AF35 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51496-1**

Capture Antibody Information

Catalog Number: 60289-2-PBS	Clone ID: 6D3D12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: U2 small nuclear RNA auxiliary factor 1
Isotype: IgG1	GenBank: BC001923	Gene ID: 7307
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag0399	

Detection Antibody Information

Catalog Number: 60289-1-PBS	Clone ID: 2B7F11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat	Full name: U2 small nuclear RNA auxiliary factor 1
Isotype: IgG2c	GenBank: BC001923	Gene ID: 7307
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0399	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51496-1 targets U2AF35 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: U2AF35 Monoclonal antibody, PBS Only (Capture) 60289-2-PBS (6D3D12). 100 µg. Concentration 1 mg/mL.

Detection antibody: U2AF35 Monoclonal antibody, PBS Only (Detector) 60289-1-PBS (2B7F11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51496-1, U2AF35 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60289-2-PBS. Detection antibody: 60289-1-PBS.



Cytometric bead array standard curve of MP51496-1, U2AF35 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60289-2-PBS. Detection antibody: 60289-1-PBS. Standard: Ag0399. Range: 3.125-400 ng/mL.