

DMXL2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51508-1**

Capture Antibody Information

Catalog Number:
66891-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B5A1
Reactivity:
human
GenBank:
BC140781
Immunogen Catalog Number:
Ag26302

Conjugate:
Unconjugated
Full name:
Dmx-like 2
Gene ID:
23312

Detection Antibody Information

Catalog Number:
66891-3-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
2F5B12
Reactivity:
human
GenBank:
BC140781
Immunogen Catalog Number:
Ag26302

Conjugate:
Unconjugated
Full name:
Dmx-like 2
Gene ID:
23312

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51508-1 targets DMXL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DMXL2 Monoclonal antibody, PBS Only (Capture) 66891-1-PBS (2B5A1). 100 µg. Concentration 1 mg/mL.

Detection antibody: DMXL2 Monoclonal antibody, PBS Only (Detector) 66891-3-PBS (2F5B12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51508-1, DMXL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66891-1-PBS. Detection antibody: 66891-3-PBS.



Cytometric bead array standard curve of MP51508-1, DMXL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66891-1-PBS. Detection antibody: 66891-3-PBS. Standard:Ag26302. Range: 1.563-200 ng/mL