

# TRIM25 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51519-1**

## Capture Antibody Information

Catalog Number:  
67314-1-PBS  
Host:  
Mouse  
Isotype:  
IgG1  
Purification Method:  
Protein A purification

Clone ID:  
1H6G1  
Reactivity:  
human  
GenBank:  
BC016924  
Immunogen Catalog Number:  
Ag29325

Conjugate:  
Unconjugated  
Full name:  
tripartite motif-containing 25  
Gene ID:  
7706

## Detection Antibody Information

Catalog Number:  
67314-2-PBS  
Host:  
Mouse  
Isotype:  
IgG1  
Purification Method:  
Protein G Magarose purification

Clone ID:  
1C5F1  
Reactivity:  
human  
GenBank:  
BC016924  
Immunogen Catalog Number:  
Ag29325

Conjugate:  
Unconjugated  
Full name:  
tripartite motif-containing 25  
Gene ID:  
7706

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51519-1 targets TRIM25 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM25 Monoclonal antibody, PBS Only (Capture) 67314-1-PBS (1H6G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM25 Monoclonal antibody, PBS Only (Detector) 67314-2-PBS (1C5F1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

## Selected Validation Data



Cytometric bead array sample test of MP51519-1, TRIM25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67314-1-PBS. Detection antibody: 67314-2-PBS.



Cytometric bead array standard curve of MP51519-1, TRIM25 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67314-1-PBS. Detection antibody: 67314-2-PBS. Standard:Ag29325. Range: 1.563-200 ng/mL