

For Research Use Only

# RBM8A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51528-1

## Capture Antibody Information

Catalog Number:  
67541-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2B11D1

Reactivity:  
human

GenBank:  
BC017088

Immunogen Catalog Number:  
Ag6800

Conjugate:  
Unconjugated

Full name:  
RNA binding motif protein 8A

Gene ID:  
9939

## Detection Antibody Information

Catalog Number:  
67541-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1A5B6

Reactivity:  
human

GenBank:  
BC017088

Immunogen Catalog Number:  
Ag6800

Conjugate:  
Unconjugated

Full name:  
RNA binding motif protein 8A

Gene ID:  
9939

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51528-1 targets RBM8A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBM8A,Y14 Monoclonal antibody, PBS Only (Capture) 67541-2-PBS (2B11D1). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBM8A,Y14 Monoclonal antibody, PBS Only (Detector) 67541-3-PBS (1A5B6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS.



Cytometric bead array standard curve of MP51528-1, RBM8A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67541-2-PBS. Detection antibody: 67541-3-PBS. Standard:Ag6800. Range: 1.563-200 ng/mL