

Ubiquilin 2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51530-1

Capture Antibody Information

Catalog Number:
60495-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3F3D8
Reactivity:
human
GenBank:
BC069237
Immunogen Catalog Number:
Ag20251

Conjugate:
Unconjugated
Full name:
ubiquilin 2
Gene ID:
29978

Detection Antibody Information

Catalog Number:
60495-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2A4A4
Reactivity:
human
GenBank:
BC069237
Immunogen Catalog Number:
Ag20251

Conjugate:
Unconjugated
Full name:
ubiquilin 2
Gene ID:
29978

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51530-1 targets Ubiquilin 2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Ubiquilin 2 Monoclonal antibody, PBS Only (Capture) 60495-2-PBS (3F3D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: Ubiquilin 2 Monoclonal antibody, PBS Only (Detector) 60495-3-PBS (2A4A4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51530-1, Ubiquilin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60495-2-PBS. Detection antibody: 60495-3-PBS.



Cytometric bead array standard curve of MP51530-1, Ubiquilin 2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60495-2-PBS. Detection antibody: 60495-3-PBS. Standard: Ag20251. Range: 0.781-100 ng/mL.