

For Research Use Only

FEZ1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51532-1

Capture Antibody Information

Catalog Number:
61033-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
6A8A8
Reactivity:
human
GenBank:
BC009545
Immunogen Catalog Number:
Ag2744

Conjugate:
Unconjugated
Full name:
fasciculation and elongation protein zeta 1 (zygin I)
Gene ID:
9638

Detection Antibody Information

Catalog Number:
61033-2-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A purification

Clone ID:
9G11A9
Reactivity:
human
GenBank:
BC009545
Immunogen Catalog Number:
Ag2744

Conjugate:
Unconjugated
Full name:
fasciculation and elongation protein zeta 1 (zygin I)
Gene ID:
9638

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51532-1 targets FEZ1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FEZ1 Monoclonal antibody, PBS Only (Capture) 61033-1-PBS (6A8A8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: FEZ1 Monoclonal antibody, PBS Only (Detector) 61033-2-PBS (9G11A9). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51532-1, FEZ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61033-1-PBS. Detection antibody: 61033-2-PBS.



Cytometric bead array standard curve of MP51532-1, FEZ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61033-1-PBS. Detection antibody: 61033-2-PBS. Standard: Ag2744. Range: 1.563-200 ng/mL.