

NR3C1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51534-1**

Capture Antibody Information

Catalog Number: 66904-2-PBS	Clone ID: 2E8B5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
Isotype: IgG1	GenBank: BC015610	Gene ID: 2908
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28431	

Detection Antibody Information

Catalog Number: 66904-3-PBS	Clone ID: 2B10B9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
Isotype: IgG1	GenBank: BC015610	Gene ID: 2908
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28431	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51534-1 targets NR3C1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Glucocorticoid receptor Monoclonal antibody, PBS Only (Capture) 66904-2-PBS (2E8B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Glucocorticoid receptor Monoclonal antibody, PBS Only (Detector) 66904-3-PBS (2B10B9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS.



Cytometric bead array standard curve of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS. Standard: Ag28431. Range: 0.781-100 ng/mL.