

ERC1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51549-2**

Capture Antibody Information

Catalog Number:
66681-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2A4B6
Reactivity:
human
GenBank:
BC132782
Immunogen Catalog Number:
Ag17665

Conjugate:
Unconjugated
Full name:
ELKS/RAB6-interacting/CAST family member 1
Gene ID:
23085

Detection Antibody Information

Catalog Number:
66681-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1G3D11
Reactivity:
human
GenBank:
BC132782
Immunogen Catalog Number:
Ag17665

Conjugate:
Unconjugated
Full name:
ELKS/RAB6-interacting/CAST family member 1
Gene ID:
23085

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51549-2 targets ERC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ERC1 Monoclonal antibody, PBS Only (Capture) 66681-3-PBS (2A4B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: ERC1 Monoclonal antibody, PBS Only (Detector) 66681-4-PBS (1G3D11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51549-2, ERCC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66681-3-PBS. Detection antibody: 66681-4-PBS.



Cytometric bead array standard curve of MP15149-2, ERC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66681-3-PBS. Detection antibody: 66681-4-PBS. Standard: Ag17665. Range: 0.781-100 ng/mL.