

SF3B3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51555-1**

Capture Antibody Information

Catalog Number:
67469-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
3G5E10
Reactivity:
human
GenBank:
BC003146
Immunogen Catalog Number:
Ag6315

Conjugate:
Unconjugated
Full name:
splicing factor 3b, subunit 3, 130kDa
Gene ID:
23450

Detection Antibody Information

Catalog Number:
67469-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3H9E3
Reactivity:
human
GenBank:
BC003146
Immunogen Catalog Number:
Ag6315

Conjugate:
Unconjugated
Full name:
splicing factor 3b, subunit 3, 130kDa
Gene ID:
23450

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51555-1 targets SF3B3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SF3B3 Monoclonal antibody, PBS Only (Capture) 67469-2-PBS (3G5E10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SF3B3 Monoclonal antibody, PBS Only (Detector) 67469-3-PBS (3H9E3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51555-1, SF3B3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67469-2-PBS. Detection antibody: 67469-3-PBS.



Cytometric bead array standard curve of MP51555-1, SF3B3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67469-2-PBS. Detection antibody: 67469-3-PBS. Standard:Ag6315. Range: 3.125-400 ng/mL.