

For Research Use Only

# TBL1X Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51581-1

## Capture Antibody Information

Catalog Number:  
66955-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1D4D12

Reactivity:  
human

GenBank:  
BC032708

Immunogen Catalog Number:  
Ag5291

Conjugate:  
Unconjugated

Full name:  
transducin (beta)-like 1X-linked

Gene ID:  
6907

## Detection Antibody Information

Catalog Number:  
66955-3-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A purification

Clone ID:  
2C2B2

Reactivity:  
human

GenBank:  
BC032708

Immunogen Catalog Number:  
Ag5291

Conjugate:  
Unconjugated

Full name:  
transducin (beta)-like 1X-linked

Gene ID:  
6907

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51581-1 targets TBL1X in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TBL1X Monoclonal antibody, PBS Only (Capture) 66955-2-PBS (1D4D12). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: TBL1X Monoclonal antibody, PBS Only (Detector) 66955-3-PBS (2C2B2). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

