

STXBP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51586-1

Capture Antibody Information

Catalog Number: 67137-2-PBS	Clone ID: 2D8B8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: syntaxin binding protein 1
Isotype: IgG3	GenBank: BC015749	Gene ID: 6812
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28605	

Detection Antibody Information

Catalog Number: 67137-3-PBS	Clone ID: 1B11A10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: syntaxin binding protein 1
Isotype: IgG1	GenBank: BC015749	Gene ID: 6812
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28605	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51586-1 targets STXBP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STXBP1 Monoclonal antibody, PBS Only (Capture) 67137-2-PBS (2D8B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: STXBP1 Monoclonal antibody, PBS Only (Detector) 67137-3-PBS (1B11A10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

