

BCAS3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51587-1**

Capture Antibody Information

| | | |
|------------------------|---------------------------|---------------------------------------|
| Catalog Number: | Clone ID: | Conjugate: |
| 67085-2-PBS | 1E1D4 | Unconjugated |
| Host: | Reactivity: | Full name: |
| Mouse | human | breast carcinoma amplified sequence 3 |
| Isotype: | GenBank: | Gene ID: |
| IgG1 | BC001250 | 54828 |
| Purification Method: | Immunogen Catalog Number: | |
| Protein G purification | Ag28545 | |

Detection Antibody Information

| | | |
|------------------------|---------------------------|---------------------------------------|
| Catalog Number: | Clone ID: | Conjugate: |
| 67085-1-PBS | 1G7B1 | Unconjugated |
| Host: | Reactivity: | Full name: |
| Mouse | human, mouse, rat, pig | breast carcinoma amplified sequence 3 |
| Isotype: | GenBank: | Gene ID: |
| IgG1 | BC001250 | 54828 |
| Purification Method: | Immunogen Catalog Number: | |
| Protein A purification | Ag28545 | |

Applications

| | | |
|-----------------------|--|--|
| Tested Applications: | Range: | Recommended Dilutions: |
| Cytometric bead array | 0.391-50 ng/mL (Cytometric Bead Array) | It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |

Product Information

MP51587-1 targets BCAS3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCAS3 Monoclonal antibody, PBS Only (Capture) 67085-2-PBS (1E1D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: BCAS3 Monoclonal antibody, PBS Only (Detector) 67085-1-PBS (1G7B1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data

