

RALBP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51591-1

Capture Antibody Information

Catalog Number: 67359-2-PBS	Clone ID: 2D1B2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: raLA binding protein 1
Isotype: IgG1	GenBank: BC013126	Gene ID: 10928
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28618	

Detection Antibody Information

Catalog Number: 67359-3-PBS	Clone ID: 2A7C9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: raLA binding protein 1
Isotype: IgG1	GenBank: BC013126	Gene ID: 10928
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28618	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51591-1 targets RALBP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RALBP1 Monoclonal antibody, PBS Only (Capture) 67359-2-PBS (2D1B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: RALBP1 Monoclonal antibody, PBS Only (Detector) 67359-3-PBS (2A7C9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.



Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51591-1, RALBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67359-2-PBS. Detection antibody: 67359-3-PBS.

Cytometric bead array standard curve of MP51591-1, RALBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67359-2-PBS. Detection antibody: 67359-3-PBS. Standard: Ag28618. Range: 3.125-400 ng/mL.