

For Research Use Only

Uteroglobin Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51593-1

Capture Antibody Information

Catalog Number:
61047-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
2C10B2
Reactivity:
human
GenBank:
BC004481
Immunogen Catalog Number:
Ag26315

Conjugate:
Unconjugated
Full name:
secretoglobin, family 1A, member 1 (uteroglobin)
Gene ID:
7356

Detection Antibody Information

Catalog Number:
61047-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1E9H6
Reactivity:
human
GenBank:
BC004481
Immunogen Catalog Number:
Ag26315

Conjugate:
Unconjugated
Full name:
secretoglobin, family 1A, member 1 (uteroglobin)
Gene ID:
7356

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51593-1 targets Uteroglobin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Uteroglobin Monoclonal antibody, PBS Only (Capture) 61047-1-PBS (2C10B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Uteroglobin Monoclonal antibody, PBS Only (Detector) 61047-2-PBS (1E9H6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51593-1, Uteroglobin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61047-1-PBS. Detection antibody: 61047-2-PBS.



Cytometric bead array standard curve of MP51593-1, Uteroglobin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61047-1-PBS. Detection antibody: 61047-2-PBS. Standard: Ag26315. Range: 0.156-20 ng/mL