

IKZF1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51595-1

Capture Antibody Information

Catalog Number:
66966-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1C1G12
Reactivity:
human
GenBank:
BC018349
Immunogen Catalog Number:
Ag28637

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 1 (Ikaros)
Gene ID:
10320

Detection Antibody Information

Catalog Number:
66966-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C9G5
Reactivity:
human
GenBank:
BC018349
Immunogen Catalog Number:
Ag28637

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 1 (Ikaros)
Gene ID:
10320

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51595-1 targets IKZF1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IKZF1 Monoclonal antibody, PBS Only (Capture) 66966-2-PBS (1C1G12). 100 µg. Concentration 1 mg/mL.

Detection antibody: IKZF1 Monoclonal antibody, PBS Only (Detector) 66966-3-PBS (2C9G5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS.



Cytometric bead array standard curve of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS. Standard: Ag28637. Range: 1.563-200 ng/mL.