

ATX3,ATXN3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51611-1

Capture Antibody Information

Catalog Number: 67057-1-PBS	Clone ID: 1F7E10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, pig	Full name: ataxin 3
Isotype: IgG2a	GenBank: BC033711	Gene ID: 4287
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6627	

Detection Antibody Information

Catalog Number: 67057-2-PBS	Clone ID: 1D11G4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ataxin 3
Isotype: IgG1	GenBank: BC033711	Gene ID: 4287
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag6627	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51611-1 targets ATX3,ATXN3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATX3,ATXN3 Monoclonal antibody, PBS Only (Capture) 67057-1-PBS (1F7E10). 100 µg. Concentration 1 mg/mL.

Detection antibody: ATX3,ATXN3 Monoclonal antibody, PBS Only (Detector) 67057-2-PBS (1D11G4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51611-1, ATX3,ATXN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67057-1-PBS. Detection antibody: 67057-2-PBS.



Cytometric bead array standard curve of MP51611-1, ATX3,ATXN3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67057-1-PBS. Detection antibody: 67057-2-PBS. Standard:Ag6627. Range: 3.125-400 ng/mL.