

For Research Use Only

WNT10B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51624-1

Capture Antibody Information

Catalog Number:
67210-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H2B3

Reactivity:
human

GenBank:
BC096353

Immunogen Catalog Number:
Ag17114

Conjugate:
Unconjugated

Full name:
wingless-type MMTV integration site family, member 10B

Gene ID:
7480

Detection Antibody Information

Catalog Number:
67210-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G2B6

Reactivity:
human

GenBank:
BC096353

Immunogen Catalog Number:
Ag17114

Conjugate:
Unconjugated

Full name:
wingless-type MMTV integration site family, member 10B

Gene ID:
7480

Applications

Tested Applications:
Cytometric bead array

Range:
3.12-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51624-1 targets WNT10B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: WNT10B Monoclonal antibody, PBS Only (Capture) 67210-2-PBS (1H2B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: WNT10B Monoclonal antibody, PBS Only (Detector) 67210-3-PBS (1G2B6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51624-1, WNT10B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67210-2-PBS. Detection antibody: 67210-3-PBS.



Cytometric bead array standard curve of MP51624-1, WNT10B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67210-2-PBS. Detection antibody: 67210-3-PBS. Standard: Ag17114. Range: 3.12-400 ng/mL.